




Introduction to FlowJo

Flow Cytometry Core Facility

Diana Ordoñez



EMBL



FlowJo software

FlowJo is an integrated environment for viewing and analyzing flow cytometric data, presented in the form of a Workspace.

The screenshot shows the FlowJo software interface. At the top is the Task Bar with icons for Open..., Print..., Save..., Revert, Save As..., Export Workspace, Data Export / Concatenate..., Find, and FCS Scan. Below the Task Bar is the ribbon menu with sections for Document, Compensation, Ab Titering, and Panel_1.2. The main workspace displays a table titled 'Group' with columns for 'Group', 'Size', and 'Role'. The table contains four entries: 'All Samples' (Size 14, Role Test), 'Compensation' (Size 0, Role Compensation), 'Ab Titering' (Size 5, Role Test), and 'Panel_1.2' (Size 9, Role Test). Below this is a hierarchical tree view of samples under 'CD8a_1,3a,100_A01.fcs'. The tree includes nodes for 'Lymphocytes', 'CD8a+', 'Median : Comp-PerCP-Cy5-5-A (CD8a)', 'CD8a-', and 'Median : Comp-PerCP-Cy5-5-A (CD8a)'. The sample table below lists five samples: CD8a_1,3a,200_A02.fcs, CD8a_1,3a,400_A03.fcs, CD8a_1,3a,800_A04.fcs, and CD8a_1,3a,1600_A05.fcs, each with collection date 27-FEB-20..., statistic values (62.5, 16.7, 11167, 83.3, 41.3), and cell counts (20000, 12502, 2082, 10420).

Group	Size	Role
{ } All Samples	14	Test
{ } Compensation	0	Compensation
► { } Ab Titering	5	Test
{ } Panel_1.2	9	Test

Name	Collecti...	Statistic	#Cells	Ab Dilution	Ab concen...
CD8a_1,3a,100_A01.fcs	27-FEB-20...		20000	1:100	100
▼ Lymphocytes		62.5	12502		
▼ CD8a+		16.7	2082		
Σ Median : Comp-PerCP-Cy5-5-A (CD8a)		11167			
▼ CD8a-		83.3	10420		
Σ Median : Comp-PerCP-Cy5-5-A (CD8a)		41.3			
► CD8a_1,3a,200_A02.fcs	27-FEB-20...		20000	1:200	200
► CD8a_1,3a,400_A03.fcs	27-FEB-20...		20000	1:400	400
► CD8a_1,3a,800_A04.fcs	27-FEB-20...		20000	1:800	800
► CD8a_1,3a,1600_A05.fcs	27-FEB-20...		20000	1:1600	1600



Task Bar

Ribbons

Group section

Sample section

Population Hierarchy





Allow the visualization of the populations defined for a sample and the relationship between gated populations

Data display

Flow cytometric data is usually display using logarithmic scale, due to the high dynamic range between negative and positive populations.



Logarithmic displays allows the visualization of subset of cells with medium to high amount of fluorescence. However, its might offer a truncated view of subsets with little, not cell-associated fluorescence or fluorescence values below zero.



Fluorescence signals below zero

Flow cytometers measure cell-associated fluorescence that can be essentially zero but not negative.

Negative – Intermediate - High



The image displayed in the graphs is no longer “pure” fluorescence , it is a corrected measurement derived from a fluorescence signal.

Negatives values are a consequence of 2 different operations related to the way the instruments collects and corrects data:



- Baseline correction
- Compensation


Bi-exponential Transformation

Transformation do not modify your actual recorded fluorescence data, only the amount of visual space that is allotted to a specific region of the data.

The biexponential scale is a hybrid of a linear and logarithmic scale. The low end of the data is displayed with a linear scale allowing negative numbers and zero to be displayed, while the high end of the data is displayed using a logarithmic scale.

Bi-exponential display helps to visualize data that is compressed against the lower x and y axis, improving it resolution.







Extra Negative Decades: Amount of negative space is shown.

The **width basis** is the amount of space shown on a linear scale, on both sides of zero. For example, if you set the width basis to -10 then the space from -10 to +10 will be shown linearly and the space above and below that will be scaled logarithmically.


Positive Decade Range (+/- buttons) - This setting allows you to control how many decades the data is displayed over.

Calculate width basis: Automatic calculation made by FlowJo.

Bi-exponential scale could be applied in the acquisition software (FacSDIVA).




no biexponential scaling



biexponential scaling (both axes)

Bi-exponential display reveals the hidden double-negative population that is not visible in the standard log scale!!!

The beginner's guide to pulse processing




Slide from Derek Davis, The Francis Crick Institute London UK.

Data Acquisition


Because fluorescent intensity is a relative measure, the position of the data in the DNA histogram can be controlled. It is convention to place the G0/G1 peak at the 50 or 100 intensity unit mark.

PMT voltage could be slightly modify during acquisition in order to place G1 peak at the same position in every sample If possible.




Flow Rate

Narrow Sample Stream:
Low Flow Rate




Wide Sample Stream:
High Flow Rate



Cell Cycle Models

FlowJo uses two model options to fit your data: Watson (Pragmatic) and the Dean-Jett-Fox models.

Comparison of the S-Phase Fit




Watson (Pragmatic) model: makes no assumptions about the shape of the S-Phase distribution.

Dean-Jet-Fox model assumes that the S-phase can be modeled by a second degree polynomial (that is convoluted with gaussian distributions of varying width throughout S-phase).

Watson (Pragmatic) Model

The Watson model assumes that only the data within the G0/G1 and G2/M peaks are normally distributed and that one of those two peaks is identifiable.




G0/G1 peaks are approximated to a Gaussian distribution

G1 mean is guess by finding the channel with the most cell in the left portion of the data.

The **standard deviation (SD)** or width of the population is approximated by finding the width of the distribution at 60% of the maximum height.

Fit improvement: A minimization process is executed over a range of -3 to 1 standard deviations about the first guess mean



G2 peak:


G2/M mean is placed at $1.75 \times$ the intensity of the G0/G1 mean.

The **standard deviation (SD)** is estimated in the same manner and a second Gaussian distribution is fit to the data using the same minimization process.

S phase:

Gaussian distributions of G0/G1 and G2/M are subtracted from the data to leave only the cells that will be considered to be S phase.



Once only S phase cells remain, a curve is fitting by an integrating function and the edges are modified in order to reflect the G0/G1 and G2/M edges.




The complete model is then the sum of the three parts.

Cell Cycle Platform

Drag the samples and start the analysis of the control sample.



Result Display




Reflects the quality of the fit.
% of cells ascribed to a given population

Gaussian distributions can be completely described by a mean (μ) and a CV.

The **Root-mean-square error (RMSE)** is a frequently used measure of the differences between values predicted by a model and the values actually observed.

Since RMS is a measure of the distance from the model to the data, a smaller RMS indicates a better fit. Usually.... RMS value of

Improving the cell cycle model



Constraint settings:

Peak location: Set G1 or G2 peak positions based on a control sample.

Mean constraint: Constraint the peak positions relative to each other.


The biology of the system tells us that G2/M cells should have twice as much DNA as G1 cells. You could define the G2 peak position to be around 2X the G1 peak or the G1 peak to be close to 0.5X the G2 peak.

CV constraint: Set the same CV values for G1 and G2 peaks.

The two populations are subject to the same illumination error derived from the analyzer.

RMS can be used to determine whether a constraint improved the fit or not.

Data display – Batch analysis



Summary table

FlowJo Tables: CC Tutorial D.wsp

Column	Population	Statistic	Parameter	Name
1 Σ	Single Cells/Cell Cycle	RMSD		
2 Σ	Single Cells/Cell Cycle	%G1		
3 Σ	Single Cells/Cell Cycle			
4 Σ	Single Cells/Cell Cycle			
5 Σ	Single Cells/Cell Cycle			
6 Σ	Single Cells/Cell Cycle			
7 Σ	Single Cells/Cell Cycle			
8 Σ	Single Cells/Cell Cycle			
9 Σ	Single Cells/Cell Cycle			
10 Σ	Single Cells/Cell Cycle			
11 Σ	Single Cells/G2			
12 Σ	Single Cells/M cells (pH3+)			

Table - Table

Ancestry Subset Statistic For	Single Cells/ Cell Cycle RMSD	Single Cells/ Cell Cycle %G1	Single Cells/ Cell Cycle %S	Single Cells/ Cell Cycle %G2	Single Cells/ Cell Cycle G1 Mean	Single Cells/ Cell Cycle G2 Mean	Single Cells/ Cell Cycle G1 CV	Single Cells/ Cell Cycle G2 CV	Single Cells/ Cell Cycle % < G1	Single Cells/ Cell Cycle % > G2
Specimen_001_HK u...	0.45	34.5	37.6	22.0	46957	91847	4.60	3.41	2.00	3.91
Specimen_001_HK n...	0.25	3.26	22.7	69.4	40601	96564	10.6	4.44	0.71	5.13
Specimen_001_HK V...	1.33	1.04	5.78	46.8	52091	99533	16.2	4.71	-0.063	44.8
Specimen_001_HK V...	0.47	0.050	44.0	44.7	24852	103323	3.32	3.32	3.00E-3	10.9
Specimen_001_HK D...	0.37	1.14	18.5	67.0	52777	105573	3.62	3.58	2.90	10.7
Specimen_001_HK 7...	0.23	33.4	36.4	20.7	55449	108193	4.18	3.11	5.40	3.85
Specimen_001_cytar...	0.13	41.7	35.2	15.2	58325	111984	4.15	3.13	3.26	4.00
Specimen_001_HK n...	0.15	2.32	27.6	54.0	52936	105872	3.45	3.45	7.46	9.09
Specimen_001_HK V...	0.24	0	8.00	34.3	56714	113405	0.24	4.48	0	56.2
Specimen_001_cytar...	0.37	50.2	33.1	15.0	61800	117188	4.86	5.11	2.08	-0.22
Specimen_001_Vincr...	0.88	1.36	22.1	64.5	52247	104543	4.35	4.28	6.17	5.47
Mean	0.44	15.4	26.5	41.2	50432	105275	5.42	3.91	2.72	14.0
SD	0.36	20.0	12.3	21.0	10184	7453	4.32	0.71	2.63	18.5

END